# uvex

## EN ISO 374

DE	Information und Gebrauchsanleitung
GB	Information and user instructions
FR	Information et notice d'utilisation
IT	Informazioni e istruzioni per l'uso
ES	Información e instrucciones de uso
PT	Informações e instruções de uso
GR	Πληροφορίες και οδηγίες χρήσης
TR	Kullanıcı talimatı ve bilgileri
DK	Information og brugsvejledning
NO	Informasjon og bruksanvisning
SE	Information och bruksanvisning
FI	Tietoja ja käyttöohje
NL	Informatie en gebruiksaanwijzing
RU	Информация и руководство по использованию
PL	Informacje i instrukcja użycia
CZ	Informace a návod k použití
SK	Informácie a návod na použitie
EE	Teave ja kasutajajuhend
LV	Informācija un lietotāja instrukcija
LT	Informacija ir naudojimo instrukcija
HU	Információ és használati utasítás
RO	Informații și instrucțiuni pentru utilizatori
BG	Информация и инструкции за употреба
SI	Informacije in navodila za uporabnike
HR	Informacije i upute za korisnika
RS	Informacije i uputstva za korisnike
JP	情報とユーザーへの指示
CN	信息和用户指南
SA	مدخلتسمل تاميلعتو تامولعمل

---

protecting people

D	Prüfinstitute	Kenn.-Nr.
GB	Approved bodies	Ref. no.
FR	Organismes CE notifiés	Réf.
IT	Organismi CE accreditati	N. rif.
ES	Entidades notificadas de la EC	Ref. No.
PT	Organismos notificados CE	Nº Ref.
GR	Κοινοποιημένοι οργανισμοί ΕΚ	Αριθ. αναφ.
TR	AT onaylı kuruluşlar	Ref. no.
DK	Bemyndigede EU-organer	Ref.nr.
NO	Teknisk kontrollorgan	Ref.nr.
SE	EG:s anmälda organ	Ref. nr
FI	EY:n ilmoitetut laitokset	Viitenumero
NL	EG-testinstituten	Referentienr.
RU	EC уполномоченные органы	Справ. no.
PL	Jednostki notyfikowane w WE	Nr ref.
CZ	Notifikované orgány ES	Ref. č.
SK	Uповедomené orgány EÚ	Ref. č.
EE	EÜ registreeritud asutused	Viitenr
LV	EK pilnvarotās iestādes	Atsauces nr.
LT	EB notifikuotoji įstaiga	Nuor. nr.
HU	EK kijelölt szerv	Száma:
RO	Organisme notificate în CE	Nr. de ref.
BG	Органи акредитирани от ЕК	Съотв. №
SI	Priglašeni organi ES	Opravlina št.
HR	Obaviještena tijela od Europske Komisije	Ref. broj:
RS	EC ovlašćena tela	Ref. br.
JP	EC 通知機関	リファレンス番号
CN	EC 公告机构	参考编号
SA	الهيئات المانحة لشهادة EC	مرجع رقم
Centexbel Gent Technologiepark 70 BE - 9052 Zwijnaarde (Gent), Belgium		0493
CTC Groupe 4. rue Hermann Frenkel 69367 Lyon Cedex 07 · France		0075
RICOTEST Via Tione 9 37010 Pastrengo (VR), Italy		0498
SATRA Technology Europe Ltd. Bracetown Business Park Clonee D15 YN2P, Ireland		2777
TÜV Rheinland LGA Products GmbH Tillystraße 2, 90431 Nürnberg, Germany		0197
BSI Group The Netherlands B.V. Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, Netherlands		2797

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## GB Information and user instructions

### Safety gloves against chemical risks

#### Pictogram and performance levels

EN ISO 374-1/Type A

J K L M N O

Specification of at least six test chemicals

EN ISO 374-1/Type B

J K L

Specification of at least three test chemicals

EN ISO 374-1/Type C

Specification of at least 1 test chemical

Safety gloves for protection against chemicals are classified into three types depending on their permeation performance:

**Type A:** The safety glove maintains permeation resistance for at least 30 minutes for at least six test chemicals.

**Type B:** The safety glove maintains permeation resistance for at least 30 minutes for at least three test chemicals.

**Type C:** The safety glove maintains permeation resistance for at least 10 minutes for at least one test chemical.

The following test chemicals are possible in accordance with EN ISO 374-1:2016+A1:2018:

Code letter	Test chemical	CAS no.
A	Methanol	67-56-1
B	Acetone	67-64-1
C	Acetonitrile	75-05-8
D	Dichloromethane	75-09-2
E	Carbon disulphide	75-15-0
F	Toluene	108-88-3
G	Diethylamine	109-89-7
H	Tetrahydrofuran	109-99-9
I	Ethyl acetate	141-78-6
J	n-Heptane	142-82-5
K	Sodium hydroxide, 40%	1310-73-2
L	Sulphuric acid, 96%	7664-93-9
M	Nitric acid, 65%	7697-37-2
N	Acetic acid, 99%	64-19-7
O	Ammonium hydroxide, 25%	1336-21-6
P	Hydrogen peroxide, 30%	7722-84-1
S	Hydrofluoric acid, 40%	7664-39-3
T	Formaldehyde, 37%	50-00-0

Sampling for penetration testing is carried out over the general test level in accordance with DIN 2859, Part 1 with AQL 1.5. (Deviating AQL values < 4.0 can be found in the brief overview)

#### Permeation / Performance Level:

Level 1 (> 10 min.)

Level 2 (> 30 min.)

Level 3 (> 60 min.)

Level 4 (> 120 min.)

Level 5 (> 240 min.)

Level 6 (> 480 min.)

## Protective gloves against dangerous chemicals and micro-organisms (EN ISO 374-5:2016)

Labelling of gloves for protection against bacteria and fungi:

EN ISO 374-5

Labelling of gloves for protection against virus, bacteria and fungi:

EN ISO 374-5

VIRUS

Resistance to penetration has been assessed under laboratory conditions and relates only to the samples tested.

**Cleaning:** The gloves can only be cleaned if no aggressive chemicals have been used. Rinse glove carefully under running, lukewarm water. Otherwise the glove is single use only.

**Warning:** Before use, the gloves should be checked for any faults or defects.

This information does not give specifications regarding the actual duration of protection at the workplace and does not distinguish between compounds and pure chemicals.

Resistance to chemicals has been assessed under laboratory conditions on samples taken solely from the palm area (with the exception of cases where the glove is 400 mm or longer – in this case, the gauntlet is also tested) and relates only to the chemicals tested. This may differ if the chemical is used in a compound.

Testing is recommended to ensure the gloves are suitable for their intended use, as the conditions at the workplace may deviate from those of the type test depending on temperature, abrasion and degradation.

If safety gloves have already been used, they may provide a lower level of resistance to hazardous chemicals due to changes to their physical properties. Due to degradation, movements, straining, friction, etc., caused by contact with chemicals, the actual period of application may be significantly reduced. In the case of aggressive chemicals, degradation may be the most significant factor to be taken into account when selecting gloves resistant to chemicals. Degradation is an adverse change to one or more properties of a material for safety gloves, caused by contact with a chemical. Signs of degradation may include scuffing, swelling, disintegration, embrittlement, discolouration, change in size, appearance, hardening and softening, etc.

All technical information relates to the delivered condition, unused and unstretched at room temperature (EN 374). Once the permeation time has expired, the glove must be disposed of.

Following contamination, a change to the performance levels specified cannot be ruled out.

Gloves should not be worn where there is a risk of being caught in moving machinery parts.

**Donning:** Grasp the glove at the bottom of the cuff and pull the glove over your hand. Align the glove on your hand. Repeat the procedure with the 2nd glove.

**Doffing:** To remove, grasp the 1st glove on the outside at the bottom of the cuff and pull the glove off over your hand. Grasp the 2nd glove on the inside of the cuff and slip it off completely. Avoid direct skin contact with the contaminated outside of the gloves.

**Application:** The duration of use depends on the degree of wear and intensity of use in the respective applications. The duration of use is influenced by the combination of several factors: Permeation due to chemicals, exposure to visible and/or UV light, exposure to high or low temperatures or changes in temperature, exposure to biological agents (pests, etc.), mechanical influences (abrasion, etc.) and contamination (splash, etc.). Discolouration, deformation and atypical material changes are indications of such factors and influences.

**Disposal:** The gloves should be disposed of properly as household waste. If they come into contact with chemicals, please consult the chemical manufacturer's waste disposal guidelines.

**Storage:** Store in a dark, dry place.

A list of potentially allergenic substances or materials that could be present is available from the manufacturer if required.

**Note:** Data regarding the legal regulations for PPE products and the Declaration of Conformity can be found in the addendum or under the product item number at the following websites:

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UVEX SAFETY Gloves GmbH & Co. KG